

BG-H

Indoor Mosquito Trap

1) trap with 5 V, 2 A power supply with 3 different attachable AC plugs (Type A, C, and G), 3 m cable, LED lights, manual in en, fr, it, es



2) 1 x BG-Sweetscent (artificial human scent, effective for 2

months)

Technical details:

Weight: 2.7 lb

Dimension: 19 x 19 x 30 cm

Trap system: 5 V

Power consumption: ~10 W



| - 1 | 1/ | n | - |
|-----|----|----|---|
| - 1 | V | IJ | |
| | J | 1- | _ |

Order No. / Description

Indoor trap for mosquitoes and other insects

10300

8G-Home

Biogents 'indoor trap BG-Home catches mosquitoes (especially Asian tiger and dengue or yellow fever mosquitoes) and other flying insects inside the house.

The trap is delivered with:

- · 1x trap body including ventilator (5 V) and EU power supply (3 m)
- · 1x BG-Homescent (odor dispenser, lasts up to 1 month)

Call Now 5503 8458

Biogents Trap System

The "Biogents Trap System" combinates the different trap types. The perfectly tuned teamwork of the Biogents Trap System guarantees the most effective success outdoors and indoors.

| Biogents Trap System Stronger together | BG-Mosquitaire + CO ₂ | BG-Mosquitaire | BG-Home | |
|---|--|---|---|-----|
| BG-Booster | BG-Booster CO. | BE-Mosquitaire | eG-Home | |
| Description | Outdoor mosquito trap | Outdoor mosquito trap | Indoor trap against mosquitoes and other flying insects | |
| Targeted mosquitoes | effective against all mosquito species incl. house and flood water mosquitoes | especially effective against tiger mosquitoes | especially effective against tiger mosquitoes | |
| Characteristics of target mosquitoes | night or day biting, wide variety in characteristics | day biting, small, black/white striped, very aggressive | day biting, small, black/white striped, very aggressive | |
| What it does | reduces biting | reduces biting | reduces biting | |
| Trap type | high performance trap | tiger mosquito trap | mosquito & insect trap | |
| Use | yard / outdoor | yard / outdoor | indoor / covered terrace | |
| Power | ✓ | ✓ | ✓ | |
| Warranty | 1 years | 1 years | 1 years | |
| In the Box | | | | |
| Trap | ✓ | ✓ | ✓ | |
| Attractant | ✓ | ✓ | ✓ | |
| CO ₂ attachment | ✓ | not needed | not needed | |
| Catch bag | ✓ | ✓ | not needed | |
| Sticky cards | not needed | not needed | not needed | |
| Features | | | | |
| Indoor | | | ✓ | |
| Yard | ✓ | ✓ | | |
| CO ₂ | ✓ | ✓ | not recommended | |
| Light | not recommend | led for outdoors | ✓ | |
| Biogents patented trap technology | ✓ | ✓ | ✓ | |
| Biogents mosquito attractant technology | ~ | V | V | |
| St. March | | | | 4.3 |



Against mosquitoes and other flying insects indoors.

The BG-Home is an excellent tool that has been specifically developed to control mosquitoes and other annoying flying insects indoors:

- · Especially attractive to the Asian tiger mosquito (Aedes albopictus) and the yellow fever (or dengue) mosquito (Aedes aegypti)
- · Does not use any toxic pesticides
- · Easy to handle


The trap's capture method is based on the same patented capture technology of all Biogents suction traps:

The BGHome imitates the manner in which odors are presented by a human body. In comparison to the outdoor suction traps, additional features are included:

- · A heated surface attracts mosquitoes
- · UV light attracts other flying insects
- · The dark appearance of the trap attracts mosquitoes

It is recommended to place the trap indoors near mosquito hiding places, such as plants or nearby where mosquitoes enter the house. Do not hide the trap, nor place it too high (not more than 1 m above the ground).

How to caught:



Biogents' mosquito trap "BG-Home" is a high-performance and eco-friendly indoor trap that targets tiger mosquitoes and other flying pest insects in the house.

Patented technology

The combination of visual cues, artificially generated air plumes, and an artificial scent that mimics human skin allows the trap to effectively imitate a human being and to attract the mosquitoes.

The artificial human scent (BG-Sweetscent) is released from a sachet that is placed inside the trap. The BG-Sweetscent is effective for two months; then it must be replaced.

A heated surface at the top of the trap provides an additional attractive stimulus. Last, but not least, the black color of the trap is also visually attractive to mosquitoes.

Approaching mosquitoes are captured with the help of a fan and drawn into the base of the trap where they dehydrate.

Other indoor insect pests can be attracted by the option to switch on UV lights.

A transparent window allows you to monitor the catch and see when it is time to empty the trap.

Without toxic ingredients

The artificial scent and the trap design are internationally patented. The trap does not require carbon dioxide, propane, butane, or other flammable gases and uses no toxic chemicals.

Product Detail:











Biogents 'Skin Odor Dispensers

Biogents 'skin odors are patented and are the result of 16. The light and permeable dispensers guarantee an years of academic research into the host-finding behavior of mosquitoes at the Zoological Department of the University of Regensburg, Germany.

optimally timed release of the scent. The odor enhances the attractiveness of the traps for mosquitoes that prefer to feed on humans, especially for tiger mosquitoes.

Order No. / Description Type

BG-Sweetscent for the BG-Mosquitaire trap +/-CO,

10536

The BG-Sweetscent is an artificial human skin odor that is formulated to last for two months. To catch tiger mosquitoes (Aedes albopictus and Aedes aegypti), the BG-Sweetscent synergistically enhances the attractiveness of the mosquito trap. For all other mosquito species, the BG-Sweetscent should be used in combination with CO₂ (BG-Mosquitaire CO₂). The BG-Sweetscent is suitable for use in Biogents'outdoor suction mosquito traps and also helps to improve the tiger mosquito catch rates of many other mosquito trap types.



Set of 3 BG-Sweetscents

10536-3

Set of 3 BG-Sweetscents for the BG-Mosquitaire trap +/- CO₂ for the whole season (six months).



BG-Homescent for the **BG-Home trap**

10301

The BG-Homescent uses the same formulation to imitate the human skin as the BG-Sweetscent. It is developed to last for one monthand to be used with the BG-Home indoor trap.



Set of 3 BG-Homescents

10301-3

Set of 3 BG-Homescentsfor the BG-Home trap (for three

months).







Replacement Parts

| Туре | Order No. / Description | |
|-------------------------------------|---|----|
| BG-Mosquitaire catch bag | 10060 The catch bag is attached to the intake funnel and catches the mosquitoes that are sucked in by the ventilator. | |
| BG-Mosquitaire funnel net | 10062 The funnel net ensures that the mosquitoes in the catch bag cannot fly or crawl out of the trap. | |
| BG-Mosquitaire gauze | 10108 The plastic gauze is made of a sturdy material that is durable, washable, and it optimizes the airflow that distributes the attractant. | |
| Y-shaped power cable, Typ II | 10541 For the BG-Mosquitaire trap. To connect multiple traps or the CCU to one power supply unit. | |
| Ventilator (fan) - IP 55, Typ II | 10055 For the BG-Mosquitaire trap. This ventilator (12 V DC, 3.6 W) is rainproof and can be submerged in salt water for hours without suffering any damage. | 23 |
| Power supply unit, Typ II | 10058 For the BG-Mosquitaire trap. Specifications: EU version*, 8 m; with 100–240 V input, 12 V DC output, and 1 A, 50–60 Hz. *Other connections are available upon request | |
| DC extension cable, Typ II | 10049 For the BG-Mosquitaire trap. To be used with Biogents power supply unit. Length: 10 m. | |
| Outdoor power supply, Typ II | 10056 For the BG-Mosquitaire trap. Specifications: EU version*, 8 m; with 100–240 V input, 12 V DC output, and 1 A, 50–60 Hz. *Other connections are available upon request | |
| BG-GAT barrier net | 10705 Barrier net with elastic band for the BG-GAT trap. It is attached to the transparent chamber of the BG-GAT. | |
| Set of 10 sticky cards | 10713 Set of 10 sticky cards for the BG-GAT trap with sticky surfaces on both sides. | |

Contact

HK Misting Limited 25/F Room S, Blk 2 Goldfield Industrial Building 144-150 Tai Lin Pai Road, Kwai Chung, NT Phone: 5169 0800 Fax: 2610 0216

E-mail: admin@hkmisting.com

Visit our webpage for manuals, reports, papers, and further information: »www.hkmisting.com

- International patents granted and pending -
- © 2020 Biogents AG & hkmisting all rights reserved -











